

SURFACE / READY TO USE DISINFECTANTS

DESCOCEPT SENSITIVE WIPES



PRESERVE

WIPES FOR ALCOHOL-BASED SURFACE DISINFECTION

1 / 10

Good material and skin compatibility
with low ethanol content

Wide range of applications
on surfaces and medical devices

Efficient disinfection and cleaning performance **with
long useful life and broad spectrum of activity**

Use safely without gloves
and save additional costs



SURFACE / READY TO USE DISINFECTANTS

DESCOCEPT SENSITIVE WIPES



PRESERVE

PRODUCT DESCRIPTION

2 / 10

Varied tasks **require**
multifaceted product solutions



DESCOCEPT SENSITIVE WIPES



DESCOCEPT SENSITIVE WIPES



DESCOCEPT SENSITIVE WIPES XL

DESCOCEPT SENSITIVE WIPES are alcoholic, ready to use wipes for disinfection and cleaning of medical devices and medical inventory as well as all kind of surfaces, especially sensitive surfaces in near patient areas.

DESCOCEPT SENSITIVE WIPES feature a patented folding technology with an optimal ratio of non-woven fleece material to impregnating solution. This ensures a quick and effective disinfection, while also being gentle on sensitive surfaces and materials.

Our **DESCOCEPT SENSITIVE WIPES XL** are convenient and time-saving. They can be used everywhere and their larger size makes them ideal for disinfecting large surfaces around the patient.

DESCOCEPT SENSITIVE WIPES are the sustainable solution made from 100 % biodegradable non-woven. The non-woven fabric, made from renewable raw materials, gives greater absorption and more efficient cleaning performance. The improved nonwoven allows optimal delivery of the impregnating solution, sufficient for an area of over 2 m² per wipe. The biodegradability of our non-woven fabrics is certified to the European standard EN 13432.

Both DESCOCEPT SENSITIVE WIPES (also with DESOTEX Bio fleece) and **DESCOCEPT SENSITIVE WIPES XL** are effective against bacteria, yeasts (including Candida auris) and enveloped viruses. **DESCOCEPT SENSITIVE WIPES** products are also suitable for reprocessing disinfectant dispensers made of anodized aluminium.

SURFACE / READY TO USE DISINFECTANTS

DESCOCEPT SENSITIVE WIPES



PRESERVE

APPLICATIONS AND NOTES

3 / 10



Range of application

Acc. to Biocidal Products Regulation (BPR)

Rapid disinfection and cleaning of alcohol-resistant medical inventory and all types of surfaces, mainly in near patient areas. Suitable for use in the food sector.

Intended purpose acc. to EU Medical Devices Regulation

For disinfection and cleaning of surfaces of non-invasive medical devices with mechanical action for professional users.

Further areas of application

In addition to the medical sector also suitable for the food sector and large canteen kitchens as well as for industry and public facilities.

Application

Wipe surfaces with DESCOCEPT SENSITIVE WIPES until completely wet and keep wet for the entire contact time. Visible soiling must be removed before disinfection.

According to the EU Medical Device Regulation, users/patients are obligated to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the EU Member State in which the user/patient is established.

Application notes

Shelf life after opening: 42 days

Composition

100 g solution contain: 45 g ethanol.

Material compatibility

Wide range of applications on surfaces and medical devices.
(see page 5 - 8)

Product status

Dual Registration (medical device/biocide)

Precautionary and hazard statements

Flammable liquid and vapour. If medical advice is needed, have product container or label at hand. Keep out of reach of children. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Keep container tightly closed. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Dispose of contents/container to approved disposal company or local collection.

For professional use only by personnel with corresponding specialist knowledge according to national directives.

Use disinfectants safely.

Always read label and product information before use.

SURFACE / READY TO USE DISINFECTANTS

DESCOCEPT SENSITIVE WIPES



PRESERVE

SPECTRUM OF EFFICACY AND CONTACT TIMES

4 / 10



SPECTRUM OF EFFICACY AND CONTACT TIMES			30 s	1 min	2 min	3 min
application recommendations for surface disinfection						
bactericidal ¹ , yeasticidal ¹	VAH ² / EN ³	with mechanical action, dirty conditions			•	
tuberculocidal	EN 14348	clean and dirty conditions		•		
mycobactericidal	EN 14348	clean and dirty conditions			•	
limited spectrum virucidal	EN 14476	clean conditions			•	
	EN 14476	dirty conditions				•
active against enveloped viruses	EN 14476	clean and dirty conditions		•		
additional results						
bactericidal	EN 16615	with mechanical action, dirty conditions			•	
	EN 13727	dirty conditions		•		
	EN 1276	dirty conditions	•			
yeasticidal	EN 16615	with mechanical action, dirty conditions			•	
	EN 13624	dirty conditions		•		
	EN 1650	dirty conditions	•			
active against adenoviruses	EN 14476	clean conditions			•	
	EN 14476	dirty conditions				•
active against noroviruses (MNV)	EN 14476	clean and dirty conditions	•			
active against vacciniaviruses	EN 14476	clean and dirty conditions		•		
active against coronaviruses	EN 14476	dirty conditions	•			

1 – including phase 2 step 1 and phase 2 step 2 tests (quantitative suspension tests and quantitative carrier tests)

2 – VAH listed 5 minutes

3 – EN 13624, EN 13727, EN 16615

SURFACE / READY TO USE DISINFECTANTS

DESCOCEPT SENSITIVE WIPES



PRESERVE

MATERIAL COMPATIBILITY

5 / 10



MATERIAL METALS	not recommended	limited recommended	recommended	APPLICATION PRODUCT EXAMPLE
Stainless steel V2A			<ul style="list-style-type: none">• Medical transport chairs• Toilet chairs	
			<ul style="list-style-type: none">• Rolling walking aids• Walking frames	
Anodized aluminium			<ul style="list-style-type: none">• Disinfectant dispensers	
Brass	•		Luer lock adapter	
	•		Stethoscope (chrome-plated brass)	
Aluminium	•			

SURFACE / READY TO USE DISINFECTANTS

DESCOCEPT SENSITIVE WIPES



PRESERVE

MATERIAL COMPATIBILITY

6 / 10



MATERIAL PLASTICS: ELASTOMERS	not recommended	limited recommended	recommended	APPLICATION PRODUCT EXAMPLE
CR (chloroprene rubber)			•	
PUR (polyurethane)			•	Medical transport chairs
EPDM (ethylene propylene diene (monomer) rubber)			•	Nursing trolleys
Silicones			•	Respiratory masks
			•	Medical keyboards and mice

SURFACE / READY TO USE DISINFECTANTS

DESCOCEPT SENSITIVE WIPES



PRESERVE

MATERIAL COMPATIBILITY

7 / 10



MATERIAL PLASTICS: THERMOPLASTICS	not recommended	limited recommended	recommended	APPLICATION PRODUCT EXAMPLE
PC (polycarbonate e.g. Makrolon)			<ul style="list-style-type: none">●	Ultrasound devices
			<ul style="list-style-type: none">●	ECG devices
			<ul style="list-style-type: none">●	EEG devices
			<ul style="list-style-type: none">●	MRT devices
			<ul style="list-style-type: none">●	CT devices
			<ul style="list-style-type: none">●	Abdominal ultrasound probes
			<ul style="list-style-type: none">●	Perfusers
			<ul style="list-style-type: none">●	HF surgery housing components
ABS (acrylonitril-butadiene-styrene)			<ul style="list-style-type: none">●	Patient monitoring monitors
			<ul style="list-style-type: none">●	Medical keyboards and mice
PE-HD (polyethylene-High Density)			<ul style="list-style-type: none">●	Storage and transport containers
PE-LD (polyethylene-Low Density)			<ul style="list-style-type: none">●	
PSU (polysulfone)			<ul style="list-style-type: none">●	Housing components
PP (polypropylene)			<ul style="list-style-type: none">●	
PUR (polyurethane)			<ul style="list-style-type: none">●	
PA (polyamide)			<ul style="list-style-type: none">●	X-ray chute
PMMA (polymethylmethacrylate)	●			Acrylic and Plexiglas, incubators

SURFACE / READY TO USE DISINFECTANTS

DESCOCEPT SENSITIVE WIPES



PRESERVE

MATERIAL COMPATIBILITY

8 / 10



MATERIAL TEXTILES				APPLICATION PRODUCT EXAMPLE
	not recommended	limited recommended	recommended	
synthetic leather (e.g. with PVC and PUR top layer)			•	Examination and treatment couches
			•	Positioning aid
			•	Gynaecological examination chair
			•	Dental chairs
			•	Operating tables
PVC (polyvinylchloride)			•	Patient transport couches
			•	floor coverings

SURFACE / READY TO USE DISINFECTANTS

DESCOCEPT SENSITIVE WIPES



PRESERVE

PURCHASING INFORMATION

9 / 10

Product	Single unit	Delivery unit	Content	Wipe size	REF
DESCOCEPT SENSITIVE WIPES	Package – DESOTEX Wipes System	9	60 Wipes	20 x 22 cm	00-323DS-To60
	Package – DESOTEX Bio Wipes	6	100 Wipes	20 x 20 cm	00-323DS-T100
DESCOCEPT SENSITIVE WIPES XL	Dispenser System – DESOTEX Wipes System	4	120 Wipes	17,5 x 36 cm	00-323DS-OSEB120

The availability of the products and container sizes depends on a completed national registration.

pH 0 1 2 3 4 5 6 7 8 9 10 11 12 13 14  CE 0482

CERTIFICATIONS



Dr. Schumacher is certified according to DIN EN 13485, DIN EN ISO 9001, DIN EN ISO 14001, BS OHSAS 18001, has a validated environment management system according to EMAS and is a member of IHO, VCI, BAH, DGSV and of the DGKH.

SURFACE / READY TO USE DISINFECTANTS

DESCOCEPT SENSITIVE WIPES



PRESERVE

PRODUCT FAMILY OVERVIEW

10 / 10



DESCOCEPT SENSITIVE WIPES



DESCOCEPT SENSITIVE WIPES



DESCOCEPT SENSITIVE WIPES XL



DESCOCEPT SENSITIVE



DESCOCEPT SENSITIVE



RECOMMENDED NON-WOVEN WIPE DISPENSER SYSTEMS



ONE SYSTEM+ PLUS



ONE SYSTEM BASIC



DESCO WIPES



ECO WIPES