

A 150

Lower basket

For the optimum loading of injector modules or inserts.



- Lower basket with two module connections
- for holding up to 2 injector modules or inserts
- Automatic closing of the connections when not in use
- Clearances: width 494, length 440 mm

EAN: 4002516854449 / Material number: 12666330 / Old Material Number: 69615001D

Product category

Washer-disinfectors, laboratory

Machine affiliation

PLW 8636 [LAB / LAB MON]

PG 8504

PG 8583

PG 8583 CD

PG 8593

PLW 8604

PLW 8683

PLW 8683 CD

PLW 8693

Product category

Basket with module connection

Lower basket

Product properties

Material

Stainless steel, Plastic

Colour

Stainless Steel, Black

Area of application

Reprocessing of pipettes/phials

Reprocessing of laboratory glassware

Reprocessing of butyrometers

Reprocessing of measuring cylinders

Dimensions and weight

External dimensions, net height in mm

154

External dimensions, net width in mm

529

External dimensions, net depth in mm

546

External dimensions, gross height in mm

175

External dimensions, gross width in mm

545

External dimensions, gross depth in mm

565

Net weight in kg

2.63

Gross weight in kg

3.29