

A 151

Lower basket

For holding standard DIN mesh trays as well as various inserts.



- Vertical clearance dependent on upper basket selected
- Loading width: 490 mm, loading depth: 498 mm

EAN: 4002516853992 / Material number: 12666340 / Old Material Number: 69615101D

Product category

| | |
|---|---|
| Fresh-water dishwashers | • |
| Washer-disinfectors, laboratory | • |
| Thermal disinfectors | • |
| Washer-disinfectors, medical | • |
| Washer-disinfectors, fire brigades and emergency services | • |

Machine affiliation

| | |
|--------------------------|---|
| PG 8058 | • |
| PG 8063 | • |
| PLW 8636 [LAB / LAB MON] | • |
| PG 8581 | • |
| PG 8562 | • |
| PG 8582 | • |
| PG 8582 CD | • |
| PG 8504 | • |
| PG 8583 | • |
| PG 8583 CD | • |
| PG 8591 | • |
| PG 8592 | • |
| PG 8593 | • |
| PLW 8604 | • |
| PLW 8683 | • |
| PLW 8683 CD | • |
| PLW 8693 | • |
| PLW 8683 [Safety] | • |
| PLW 8693 [Safety] | • |
| PWA 8672 | • |
| PWD 8682 | • |
| PWD 8682 CD | • |
| PWD 8692 | • |

Product category

| | |
|------------------------|---|
| Carrier/basket carrier | • |
| Lower basket | • |

Product properties

| | |
|----------|--------------------------|
| Material | Stainless steel, Plastic |
| Colour | Stainless Steel, Black |

Area of application

| | |
|--------------------------------------|---|
| Reprocessing of laboratory glassware | • |
| Reprocessing of GYN instruments | • |
| Reprocessing of ENT instruments | • |
| Reprocessing of dental instruments | • |
| Reprocessing of theatre shoes | • |
| Reprocessing of basins | • |
| Reprocessing of kidney dishes | • |
| Reprocessing of baby bottles | • |
| Reprocessing of surgical instruments | • |
| Reprocessing of mesh trays | • |
| Reprocessing of different utensils | • |

A 151

Lower basket

For holding standard DIN mesh trays as well as various inserts.



EAN: 4002516853992 / Material number: 12666340 / Old Material Number: 69615101D

| | |
|--|------|
| Reprocessing of breathing air regulators | • |
| Reprocessing of breathing masks | • |
| Reprocessing of compressed air unit carrier frames | • |
| Capacity | |
| DIN mesh trays [number] | 2 |
| SPRI mesh trays [number] | 1 |
| Dimensions and weight | |
| External dimensions, net height in mm | 88 |
| External dimensions, net width in mm | 529 |
| External dimensions, net depth in mm | 522 |
| External dimensions, gross height in mm | 120 |
| External dimensions, gross width in mm | 545 |
| External dimensions, gross depth in mm | 535 |
| Net weight in kg | 1.84 |
| Gross weight in kg | 2.40 |