

A 202

Mobile unit

for the optimum loading of 4 DIN mesh trays.



- Suitable for 4 DIN mesh trays on 2 levels
- Insert frames: Removable for the reprocessing of tall wash items
- 6 additional connection options via Luer lock connection
- With integrated spray arm
- Loading volumes for the different levels differ from each other

EAN: 4002516854043 / Material number: 12667130 / Old Material Number: 69620201D

Product category	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Mobile unit for inserts	•
Mobile unit for mesh trays	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of surgical instruments	•
Capacity	
DIN mesh trays [number]	4
Dimensions and weight	
External dimensions, net height in mm	223
External dimensions, net width in mm	529
External dimensions, net depth in mm	542
External dimensions, gross height in mm	270
External dimensions, gross width in mm	570
External dimensions, gross depth in mm	550
Net weight in kg	6.04
Gross weight in kg	6.76