


A 302/3

Injector module

for optimum loading of laboratory glassware – can be equipped with 32 nozzles.



- Suitable for lab. glassware incl. Erlenmeyer and round flasks.
- Equipped with 16 A 844 nozzles and 16 A 845 nozzles
- Part of the EasyLoad system
- Nozzle spacing 60 mm across the width, 57 mm across the depth

EAN: 4002516188988 / Material number: 11116960 / Old Material Number: 69630203D

Product category

Washer-disinfectors, laboratory

•

Large-capacity washer-disinfectors, laboratory

•

Machine affiliation

PLW 8636 [LAB / LAB MON]

•

PG 8504

•

PG 8583

•

PG 8583 CD

•

PG 8593

•

PLW 8604

•

PLW 8683

•

PLW 8683 CD

•

PLW 8693

•

PLW 7111

•

Product category

Insert/module for narrow-necked glassware

•

Insert for various flasks

•

Product properties

Material

Stainless steel, Plastic

Colour

Stainless Steel, Black

Part of the EasyLoad system

•

Area of application

Reprocessing of laboratory glassware

•

Capacity

Narrow-necked glassware [number]

32

Dimensions and weight

External dimensions, net height in mm

143

External dimensions, net width in mm

238

External dimensions, net depth in mm

479

External dimensions, gross height in mm

135

External dimensions, gross width in mm

290

External dimensions, gross depth in mm

590

Net weight in kg

1.49

Gross weight in kg

2.67

Companion pack

A 844, injector nozzle Ø 2.5 mm, length 80 mm [number]

16

A 845, injector nozzle Ø 2.5 mm, length 125 mm [number]

16