


A 303
Module
for the optimum loading of up to 98 pipettes.



- Suitable for 98 measuring and full pipettes, for example
- Height of holding frame 150 mm
- 1 unit can be inserted into each of the baskets A 150 and APLW 106
- Vertical clearance without upper basket: 450 mm (with A 150 lower basket)

EAN: 4002516854746 / Material number: 12669220 / Old Material Number: 69630301D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8636 [LAB / LAB MON]	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
Product category	
Insert/module for pipettes	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of pipettes/phials	•
Capacity	
Pipettes [number]	98
Dimensions and weight	
External dimensions, net height in mm	185
External dimensions, net width in mm	222
External dimensions, net depth in mm	471
External dimensions, gross height in mm	220
External dimensions, gross width in mm	260
External dimensions, gross depth in mm	490
Net weight in kg	2.50
Gross weight in kg	3.08