


A 304 Module

For the optimum loading of up to 98 vials.



- Suitable for 98 centrifuge tubes, phials, test tubes, for example
- The loading heights differ depending on the positioning
- 1 unit can be inserted into each of the baskets A 100, A 150 and APLW 106

EAN: 4002516854982 / Material number: 12669230 /
Old Material Number: 69630401D

Product category

Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•

Machine affiliation

PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•

Product category

Insert/module for test tubes	•
------------------------------	---

Product properties

Material	Stainless steel, Plastic
Colour	Stainless Steel, Black

Area of application

Reprocessing of pipettes/phials	•
Reprocessing of test tubes	•

Capacity

Phials [number]	98
-----------------	----

Dimensions and weight

External dimensions, net height in mm	130
External dimensions, net width in mm	222
External dimensions, net depth in mm	471
External dimensions, gross height in mm	150
External dimensions, gross width in mm	230
External dimensions, gross depth in mm	490
Net weight in kg	2.30
Gross weight in kg	3.15