


A 11/1

Insert

For the gentle cleaning of delicate items by interrupting the spray jet.

EAN: 4002516854913 / Material number: 12666290 / Old Material Number: 69501104D



- With perforations 7 x 7 mm
- For use in upper and lower basket
- Recommended for use in PG 8536

Product category	
Washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
Product category	
Retaining frame	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	1
External dimensions, net width in mm	433
External dimensions, net depth in mm	433
External dimensions, gross height in mm	20
External dimensions, gross width in mm	470
External dimensions, gross depth in mm	470
Net weight in kg	0.70
Gross weight in kg	0.90