


E 105/1

Insert

For the optimum loading of test tubes up to 12 x 165 mm.



- Suitable for about 200 Test tubes (up to 12 x 165 mm)
- Capacity including lid A 13
- The mesh width is 9 x 9 mm

EAN: 4002516853978 / Material number: 12667380 / Old Material Number: 69510503D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
Product category	
Insert/module for test tubes	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of test tubes	•
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	200
Dimensions and weight	
External dimensions, net height in mm	212
External dimensions, net width in mm	200
External dimensions, net depth in mm	210
External dimensions, gross height in mm	235
External dimensions, gross width in mm	240
External dimensions, gross depth in mm	250
Net weight in kg	0.88
Gross weight in kg	1.15
Companion pack	
A 13 lid [number]	1