


E 144

Insert

For the optimum loading of up to 18 beakers.



- Suitable for 18 glass beakers
- equipped with 18 x 3 stay bars
- For use in upper and lower basket

EAN: 4002511115996 / Material number: 03808710 / Old Material Number: 69514401D


Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
SlimLine laboratory washers	•
Cleaning machines for industrial parts cleaning	•
Large-chamber laboratory washers for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
Product category	
Insert/module for beakers	•
Product properties	
Material	Stainless steel, Wire, rilsan-coated
Colour	Stainless Steel, White
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Beakers [quantity]	18
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	131
External dimensions, net width in mm	200
External dimensions, net depth in mm	445

E 144

Insert

For the optimum loading of up to 18 beakers.

EAN: 4002511115996 / Material number: 03808710 / Old Material Number: 69514401D



External dimensions, gross height in mm	180
External dimensions, gross width in mm	250
External dimensions, gross depth in mm	490
Net weight in kg	1.37
Gross weight in kg	1.38