


A 103

Upper basket

For the optimum loading of inserts.



- Vertical clearance: 95 mm, width 505 mm, depth 472 mm
- Particularly suitable for the reprocessing of DIN mesh trays
- Also suitable for further applications with low height items
- Can be used in combination with the A 202 mobile unit
- With integrated spray arm

EAN: 4002516855040 / Material number: 12664750 / Old Material Number: 69610302D

Product category	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PLW 8636 [LAB / LAB MON]	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 8683 [Safety]	•
PLW 8693 [Safety]	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Carrier/basket carrier	•
Upper basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Maximum vertical clearance in mm	95
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of theatre shoes	•
Reprocessing of surgical instruments	•
Reprocessing of mesh trays	•
Reprocessing of different utensils	•
Capacity	
DIN mesh trays [number]	2
SPRI mesh trays [number]	1
Dimensions and weight	
External dimensions, net height in mm	133
External dimensions, net width in mm	528
External dimensions, net depth in mm	528
External dimensions, gross height in mm	170
External dimensions, gross width in mm	545
External dimensions, gross depth in mm	560
Net weight in kg	4.18
Gross weight in kg	4.64