

A 208

Mobile unit

For the optimum loading of 8 mesh trays.



- High capacity: simultaneous reprocessing of 8 DIN mesh trays
- 4 additional connections for hollow instruments
- Removable level increases vertical clearance if required

EAN: 4002516854494 / Material number: 12668770 / Old Material Number: 69620801D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Mobile unit for mesh trays	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of mesh trays	•
Capacity	
DIN mesh trays [number]	8
Lumened instruments [number]	4
Dimensions and weight	
External dimensions, net height in mm	397
External dimensions, net width in mm	529
External dimensions, net depth in mm	546
External dimensions, gross height in mm	570
External dimensions, gross width in mm	480
External dimensions, gross depth in mm	540
Net weight in kg	9.83
Gross weight in kg	11.02