

A 313

Module

For reprocessing large-volume laboratory glassware.



- Laboratory flasks, up to 10 l
- Round flasks, up to 6 l
- Erlenmeyer flasks, up to 5 l
- Measuring flasks, up to 5 l
- Can be equipped with nozzles & holders based on individual needs
- Equipped with 5 injector positions for 2 or 3 load items

EAN: 4002516087854 / Material number: 11055970 / Old Material Number: 69631300D

Product category

Washer-disinfectors, laboratory

•

Large-capacity washer-disinfectors, laboratory

•

Machine affiliation

PLW 8636 [LAB / LAB MON]

•

PG 8504

•

PG 8583

•

PG 8583 CD

•

PG 8593

•

PLW 8604

•

PLW 8683

•

PLW 8683 CD

•

PLW 8693

•

PLW 7111

•

Product category

Insert/module for narrow-necked glassware

•

Module for large-volume laboratory glassware

•

Product properties

Material

Stainless steel, Plastic

Colour

Stainless Steel, Black

Area of application

Reprocessing of laboratory glassware

•

Dimensions and weight

External dimensions, net height in mm

221

External dimensions, net width in mm

234

External dimensions, net depth in mm

495

External dimensions, gross height in mm

260

External dimensions, gross width in mm

270

External dimensions, gross depth in mm

505

Net weight in kg

2.82

Gross weight in kg

3.29