


A 500
Mobile unit
with 1 level for cleaning laboratory glassware and utensils.



- For holding inserts and mesh trays
- Extendable to 2 levels with A 603

EAN: 4002516855026 / Material number: 12684430 /
Old Material Number: 69650002D

Product category

| | |
|---|---|
| Large-capacity washer-disinfectors, laboratory | • |
| Large-chamber washer-disinfectors, fire brigades and emergency services | • |

Machine affiliation

| | |
|-------------------|---|
| PLW 8615 | • |
| PLW 8615 [Safety] | • |
| PLW 8616 | • |
| PLW 8616 [Safety] | • |
| PLW 8617 | • |
| PLW 8617 [Safety] | • |

Product category

| | |
|-------------------------|---|
| Mobile unit for inserts | • |
|-------------------------|---|

Product properties

| | |
|----------|--------------------------|
| Material | Stainless steel, Plastic |
| Colour | Stainless Steel, Black |

Area of application

| | |
|--------------------------------------|---|
| Reprocessing of laboratory glassware | • |
|--------------------------------------|---|

Capacity

| | |
|--------------------------|---|
| DIN mesh trays [number] | 3 |
| SPRI mesh trays [number] | 2 |

Dimensions and weight

| | |
|---|------|
| External dimensions, net height in mm | 153 |
| External dimensions, net width in mm | 650 |
| External dimensions, net depth in mm | 790 |
| External dimensions, gross height in mm | 190 |
| External dimensions, gross width in mm | 655 |
| External dimensions, gross depth in mm | 810 |
| Net weight in kg | 6.02 |
| Gross weight in kg | 7.36 |