

A 501
Mobile unit
with 4 levels for cleaning laboratory glassware and utensils.



- For holding inserts and mesh trays
- Extendable to 5 or 6 levels with A 601 racks

EAN: 4002516855361 / Material number: 12684440 / Old Material Number: 69650102D

Product category

Large-capacity washer-disinfectors, laboratory

Machine affiliation

PLW 8615

PLW 8616

PLW 8617

Product category

Mobile unit for inserts

Product properties

Material

Stainless steel, Plastic

Colour

Stainless Steel, Black

Area of application

Reprocessing of laboratory glassware

Capacity

DIN mesh trays [number]

18

SPRI mesh trays [number]

12

Dimensions and weight

External dimensions, net height in mm

659

External dimensions, net width in mm

650

External dimensions, net depth in mm

710

External dimensions, gross height in mm

710

External dimensions, gross width in mm

800

External dimensions, gross depth in mm

860

Net weight in kg

21.00

Gross weight in kg

32.80