

A 503

Mobile unit

no levels, for cleaning laboratory glassware/utensils.



- For loading modules and inserts
- Depending on the height, up to 6 modules can be used
- Upper rail and upper strut can be removed

EAN: 4002516854784 / Material number: 12685380 /
Old Material Number: 69650302D

Product category

Large-capacity washer-disinfectors, laboratory

Machine affiliation

PLW 8615

PLW 8616

PLW 8617

Product category

Mobile unit with module connection

Product properties

Material

Stainless steel, Plastic

Colour

Stainless Steel, Black

Area of application

Reprocessing of laboratory glassware

Dimensions and weight

External dimensions, net height in mm

433

External dimensions, net width in mm

790

External dimensions, net depth in mm

650

External dimensions, gross height in mm

790

External dimensions, gross width in mm

715

External dimensions, gross depth in mm

860

Net weight in kg

11.40

Gross weight in kg

22.20