

A 505
Mobile unit
with 3 levels for cleaning laboratory glassware and utensils.



- For holding inserts and mesh trays

EAN: 4002516854678 / Material number: 12685400 / Old Material Number: 69650501D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Mobile unit for inserts	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of mesh trays	•
Reprocessing of different utensils	•
Capacity	
DIN mesh trays [number]	9
SPRI mesh trays [number]	6
Dimensions and weight	
External dimensions, net height in mm	543
External dimensions, net width in mm	650
External dimensions, net depth in mm	790
External dimensions, gross height in mm	790
External dimensions, gross width in mm	715
External dimensions, gross depth in mm	860
Net weight in kg	15.60
Gross weight in kg	26.40