

A 620

Injector module

for the optimum reprocessing of laboratory glassware, 10 injector positions.



- Can be equipped with up to 10 nozzles based on individual needs
- For use in A 503 mobile unit
- Part of the EasyLoad system

EAN: 4002516853909 / Material number: 12678100 / Old Material Number: 69662001D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Insert/module for narrow-necked glassware	•
Insert for various flasks	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Part of the EasyLoad system	•
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	113
External dimensions, net width in mm	133
External dimensions, net depth in mm	608
External dimensions, gross height in mm	145
External dimensions, gross width in mm	150
External dimensions, gross depth in mm	620
Net weight in kg	1.59
Gross weight in kg	1.82