

A 623  
Module  
for the optimum loading of measuring cylinders.



- Equipped with 4 nozzles
- For measuring cylinders with 1–2 l capacity
- Pivoting supports for ease of use

EAN: 4002516853893 / Material number: 12678510 / Old Material Number: 69662301D

## Product category

Large-capacity washer-disinfectors, laboratory

## Machine affiliation

PLW 8615

PLW 8616

PLW 8617

## Product category

Insert/module for measuring cylinder

## Product properties

Material

Stainless steel, Plastic

Colour

Stainless Steel, Black

## Area of application

Reprocessing of laboratory glassware

Reprocessing of measuring cylinders

## Capacity

Measuring cylinders [number]

4

## Dimensions and weight

External dimensions, net height in mm

437

External dimensions, net width in mm

255

External dimensions, net depth in mm

608

External dimensions, gross height in mm

269

External dimensions, gross width in mm

400

External dimensions, gross depth in mm

615

Net weight in kg

2.42

Gross weight in kg

3.42