

A 840/1

Injector nozzle

with plastic supports, nozzle diameter 6 mm, length 130 mm, 5 items.



- Suitable for mobile injector units & baskets screw connection
- Ideal for load items with neck diameters of between 12 mm and 85 mm
- Diameter of cap at top: 11 mm
- Part of the EasyLoad system

EAN: 4002516188803 / Material number: 11112600 / Old Material Number: 69684001D

Product category

Washer-disinfectors, laboratory

•

Large-capacity washer-disinfectors, laboratory

•

Machine affiliation

PLW 8636 [LAB / LAB MON]

•

PG 8504

•

PG 8583

•

PG 8583 CD

•

PG 8593

•

PLW 8604

•

PLW 8615

•

PLW 8616

•

PLW 8617

•

PLW 8683

•

PLW 8683 CD

•

PLW 8693

•

PLW 7111

•

Product category

Injector nozzle

•

Product properties

Material

Stainless steel, Plastic

Colour

Stainless Steel, Black

Part of the EasyLoad system

•

Area of application

Reprocessing of laboratory glassware

•

Dimensions and weight

External dimensions, net height in mm

130

External dimensions, net width in mm

92

External dimensions, net depth in mm

86

External dimensions, gross height in mm

107

External dimensions, gross width in mm

130

External dimensions, gross depth in mm

195

Net weight in kg

0.13

Gross weight in kg

0.24