

A 845

Injector nozzle

with plastic supports, nozzle diameter 2.5 mm, length 125 mm, 1 item.



- Suitable for mobile injector units & baskets screw connection
- Ideal for load items with neck diameters of between 6 mm and 55 mm
- Diameter of cap at top: 5 mm
- Part of the EasyLoad system

EAN: 4002516188797 / Material number: 11056670 / Old Material Number: 69684500D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Part of the EasyLoad system	•
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	125
External dimensions, net width in mm	60
External dimensions, net depth in mm	55
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Net weight in kg	0.01
Gross weight in kg	0.01