

A 846

Injector nozzle

with support, for cleaning laboratory flasks and round flasks.



- For laboratory bottles, max. 10 l, and round flasks, max. 6 l
- Nozzle length: 127 mm, nozzle length with holder: 158 mm
- Nozzle Ø: 10 mm
- Suitable for A 313 and A 612 modules

EAN: 4002516855088 / Material number: 12665100 / Old Material Number: 69684601D

Product category

Washer-disinfectors, laboratory

•

Large-capacity washer-disinfectors, laboratory

•

Machine affiliation

PLW 8636 [LAB / LAB MON]

•

PG 8504

•

PG 8583

•

PG 8583 CD

•

PG 8593

•

PLW 8604

•

PLW 8683

•

PLW 8683 CD

•

PLW 8693

•

PLW 7111

•

Product category

Injector nozzle

•

Product properties

Material

Stainless steel, Rubber, Plastic

Colour

Stainless Steel, Black, Grey

Area of application

Reprocessing of laboratory glassware

•

Dimensions and weight

External dimensions, net height in mm

158

External dimensions, net width in mm

110

External dimensions, net depth in mm

110

External dimensions, gross height in mm

160

External dimensions, gross width in mm

145

External dimensions, gross depth in mm

265

Net weight in kg

0.23

Gross weight in kg

0.38