

A 321
Insert
for loading different utensils.



- Suitable for beakers, funnels, breast milk pumps, for example
- With flexible handles

EAN: 4002516237099 / Material number: 11309380 / Old Material Number: 69632100D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 7111	•
Product category	
Insert/module for beakers	•
Insert/module for wide-necked glassware	•
Other inserts	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	81
External dimensions, net width in mm	240
External dimensions, net depth in mm	545
External dimensions, gross height in mm	90
External dimensions, gross width in mm	255
External dimensions, gross depth in mm	585
Net weight in kg	0.76
Gross weight in kg	0.99