

## APFD 200

### Use

for the optimum loading of air regulators.



- Suitable for holding 6 air regulators securely
- Accessories: 2 small-part baskets with 3 coloured compartments
- Breathing air connection and compressed air indicator
- Connection via quick couplings (compatible with Parker Rectus 25)
- Max. load item clearance 275 mm, diameter 12–95 mm

EAN: 4002516855330 / Material number: 12615340 / Old Material Number: 67020002D

#### Product category

Fresh-water dishwashers

•

Washer-disinfectors, fire brigades and emergency services

•

#### Machine affiliation

PG 8063

•

PLW 8683 [Safety]

•

PLW 8693 [Safety]

•

#### Product category

Insert for air regulators

•

#### Product properties

Material

Stainless steel

Colour

Stainless Steel

Compressed air indicator

•

#### Area of application

Reprocessing of breathing air regulators

•

#### Capacity

Air regulators [number]

6

#### Dimensions and weight

External dimensions, net height in mm

153

External dimensions, net width in mm

484

External dimensions, net depth in mm

452

External dimensions, gross height in mm

19

External dimensions, gross width in mm

500

External dimensions, gross depth in mm

480

Net weight in kg

7.05

Gross weight in kg

7.61