


APLW 106

Basket

for optimum loading of injector modules or inserts.



- For use on levels 1 to 4
- Basket with two docking connections
- For loading up to 2 injector modules or inserts
- Automatic closing of the connections when not in use
- Clearances dependent on which baskets are used and their position

EAN: 4002516524717 / Material number: 11854610 / Old Material Number: 69610600D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 7111	•
Product category	
Basket with module connection	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Part of the EasyLoad system	•
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	151
External dimensions, net width in mm	559
External dimensions, net depth in mm	621
External dimensions, gross height in mm	175
External dimensions, gross width in mm	570
External dimensions, gross depth in mm	630
Net weight in kg	3.05
Gross weight in kg	3.75