


## APLW 107

### Basket

height adjustable, for optimum loading of inserts.



- For use on levels 2 to 4
- Open front
- For optimum loading of inserts
- With integrated spray arm
- Clearances dependent on which baskets are used and their position

EAN: 4002516525998 / Material number: 11854290 / Old Material Number: 69610700D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Large-chamber washer-disinfectors, fire brigades and emergency services	•
Machine affiliation	
PLW 7111	•
PLW 7111 [Safety]	•
Product category	
Carrier/basket carrier	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Height-adjustable	•
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	151
External dimensions, net width in mm	559
External dimensions, net depth in mm	594
External dimensions, gross height in mm	180
External dimensions, gross width in mm	575
External dimensions, gross depth in mm	650
Net weight in kg	3.39
Gross weight in kg	5.01