


APLW 511

Mobile unit

no levels, for cleaning laboratory glassware/utensils.



- Depending on the height, up to 9 modules can be used
- Some rails and struts can be removed

EAN: 4002516525974 / Material number: 11852850 / Old Material Number: 69650300D

Product category

- Large-capacity washer-disinfectors, laboratory
- SlimLine laboratory washers for industrial parts cleaning

Machine affiliation

- PLW 8615
- PLW 8616
- PLW 8617

Product category

- Mobile unit with module connection

Product properties

- Material: Stainless steel, Plastic
- Colour: Stainless Steel, Black
- Part of the EasyLoad system

Area of application

- Reprocessing of laboratory glassware

Dimensions and weight

- External dimensions, net height in mm: 534
- External dimensions, net width in mm: 636
- External dimensions, net depth in mm: 790
- External dimensions, gross height in mm: 800
- External dimensions, gross width in mm: 720
- External dimensions, gross depth in mm: 860
- Net weight in kg: 20.16
- Gross weight in kg: 31.20