

DOS K 85/1 flex

Dispenser module for 5 litre canister

with level indicator for the automatic dispensing of liquid media.

EAN: 4002516365297 / Material number: 11498920 / Old Material Number: 69747053D



- Dispenser monitoring – high process security as per EN ISO 15883
- For dispensing liquid cleaning agent/neutralising agent
- 214 mm mm siphon for 5 l media canisters with flexible cap
- Siphon conversion kit for 10–30 l media canisters available

Product category	
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8583	•
PG 8591	•
PG 8592	•
PG 8593	•
Product category	
Dispenser system	•
Product properties	
Production rate in ml/60 s	55
Length of suction lance in mm	214
Standard electrical connection	
Number of phases	1
Neutral conductor	N
AC/DC	AC
Voltage in V	200-230
Frequency in Hz	50
Total rated load in kW	0.01
Dimensions and weight	
External dimensions, gross height in mm	160
External dimensions, gross width in mm	198
External dimensions, gross depth in mm	478
Net weight in kg	1.52
Gross weight in kg	1.74
Standards, test certificates and symbols	
CE	•