


E 109

Insert

for the optimum loading of beakers up to 250 ml.



- Suitable for 21 glass beakers, max. 250 ml
- equipped with 21 x 3 stay bars
- For use in upper and lower basket

EAN: 4002511115729 / Material number: 03808360 / Old Material Number: 69510901D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
Product category	
Insert/module for beakers	•
Product properties	
Material	Stainless steel, Wire, rilsan-coated
Colour	Stainless Steel, White
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Beakers [quantity]	21
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	155
External dimensions, net width in mm	230
External dimensions, net depth in mm	460
External dimensions, gross height in mm	185
External dimensions, gross width in mm	230
External dimensions, gross depth in mm	490
Net weight in kg	1.21
Gross weight in kg	1.54