

E 142
Insert
For holding different instruments.



- 1 mm wire, 5 mm mesh size, 5 mm frames
- equipped with 2 hinged handles
- Maximum load 10 kg

EAN: 4002511115972 / Material number: 03830300 / Old Material Number: 69514201D

Product category	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7824	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Mesh tray	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of surgical instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	55
External dimensions, net width in mm	255
External dimensions, net depth in mm	480
External dimensions, gross height in mm	75
External dimensions, gross width in mm	280
External dimensions, gross depth in mm	515

E 142

Insert

For holding different instruments.

EAN: 4002511115972 / Material number: 03830300 / Old Material Number: 69514201D

Net weight in kg	0.70
------------------	------

Gross weight in kg	0.86
--------------------	------

