

E 373

Insert

For the optimum loading of ENT instruments.



- Mesh width on the base 3 mm, sides 1.7 mm, lid 3 mm
- equipped with 28 vertical stay bars
- equipped with 2 hinged handles
- For use in upper and lower basket

EAN: 400251116405 / Material number: 03809790 / Old Material Number: 69537301D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8534	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Insert for ENT instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of ENT instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	55
External dimensions, net width in mm	150
External dimensions, net depth in mm	225
External dimensions, gross height in mm	80
External dimensions, gross width in mm	170
External dimensions, gross depth in mm	240
Net weight in kg	0.66
Gross weight in kg	0.77